



Ni Smart Magarose Beads

Index

1. Product Description.....	1
2. Purification Procedure.....	1
3. Additional Information.....	2
4. Related Products.....	2

1. Product Description

Ni Smart Magarose Beads is a new kind of IMAC packing, chelating nickel ions (Ni^{2+}) solid. **Ni Smart Magarose Beads** can be used to purify 6xHis-tagged proteins expressed in series of expression vectors, such as *E.coli.*, yeast, insect cells and mammalian cells. This form is very stable octahedral structure of nickel ions in the center, which can protect the nickel ions from attack of the competitive small molecule. The structure of Ni-SMART is compatible with a certain concentration of reducing agents, denaturing agents, detergents and other additives.

Table 1. Characteristics of **Ni Smart Magarose Beads**

Item	Description
Matrix Spherical	Magnetic agarose
Static Binding Capacity	> 10mg 6×His-tagged protein/ml medium
Particle size (μm)	30-100um
Beads concentration	20%(V/V) slurry
Storage Solution	1×PBS containing 20% ethanol
Storage Temperature	2-8 °C

Table 2. Chemical compatibilities for **Ni Smart Magarose Beads**

Reagent	Stability
0.01M NaOH	1week
10mM EDTA, 1M NaOH, 5mM DTT, 5mM TCEP, 20mM β-mercaptoethanol, 6MGua-HCl	24hour
500mM imidazole, 100mM EDTA	2hour
30% isopropanol	20min

2. Purification Procedure

2.1 Buffer Preparation

The basic principle of the following recommended buffer and other buffer is low concentration of imidazole in Lysis and wash buffer and high in elution buffer. Water and chemicals used for buffer preparation should be of high purity. It is recommended to filter the buffers by passing them through a 0.22μm or 0.45 μm filter before use.

Lysis Buffer: 50 mM NaH_2PO_4 , 0.5M NaCl, pH7.4

Wash Buffer: 50 mM NaH_2PO_4 , 0.5M NaCl, 0-5mM imidazole, pH7.4

Elution Buffer: 50 mM NaH_2PO_4 , 0.5M NaCl, 250mM imidazole, pH7.4

2.2 Sample Preparation

- 1) Centrifuge the homogenized lysate at 7,000rpm(7,500×g) for 20min at 4 °C to clarify sample. Save supernatant.
- 2) The sample contains a tolerable range of reagents and does not require dialysis or concentration. It can be sampled directly.

2.3 Preparation of the Magnetic Beads

The protocol uses 100 μl **Ni SMART Magarose Beads**, but this may be scaled up or down as required.

- 1) Completely resuspend the beads by shaking or vortexing the vial.
- 2) Transfer 500 μl **Ni Smart Magarose Beads** (20% v/v) into a clean tube.
- 3) Place the tube on a magnetic separation rack to collect the beads. Remove and discard the supernatant.
- 4) Add 1 ml Lysis Buffer to the tube and invert the tube several times to mix. Use the magnetic separation rack to collect the beads and discard the supernatant. Repeat this step twice.

2.4 Protein Purification

- 1) Transfer the 6X His-tagged Protein lysate to **Ni Smart Magarose Beads**. Incubate the tube at room temperature with mixing (on a shaker or rotator) for 10 -30 minutes. Place the tube on a magnetic separation rack to collect the beads. Remove and discard the supernatant. If necessary, keep the supernatant for analysis.
- 2) Add 1ml Wash Buffer to the tube and mix well, use the magnetic separation rack to collect the beads and discard the supernatant. Repeat the wash step three more times.
- 3) Add 200-300µl Elution Buffer, mix well and incubation for 5-10min. Use the magnetic separation rack to collect the beads and transfer the supernatant into a new tube. Repeat this step one more time.

2.5 Analysis

Identify the fractions containing the His-tagged protein. Use UV absorbance, SDS-PAGE, or western blot.

3. Additional Information

- 1) Please read the product instruction carefully before using the product.
- 2) In the process of magnetic beads preservation, operations such as freezing, drying and high-speed centrifugation should be avoided, otherwise the structure of magnetic beads will be damaged and the binding capacity of proteins will be seriously affected.
- 3) Before using the magnetic beads, please oscillate gently and fully to keep the beads in a uniform suspension state.
- 4) The beads can be reused to purify the same protein; when different proteins are purified, it is recommended to use new magnetic beads to avoid cross-contamination.

4. Related Products

Product	Cat. No.	Size
Ni IDA Magarose Beads	SM001001	1 ml
	SM001005	5 ml
	SM001025	25 ml
	SM001100	100 ml
	SM00101L	1 L
Ni NTA Magarose Beads	SM008001	1 ml
	SM008005	5 ml
	SM008025	25 ml
	SM008100	100 ml
	SM00801L	1 L
Ni Smart Magarose Beads	SM025001	1 ml
	SM025005	5 ml
	SM025025	25 ml
	SM025100	100 ml
	SM0501L	1 L